



## Aruba Health Distributors

*Aruba Health Distributors*

## DermoTherap® Abrader

[Call for Pricing](#)

### DermoTherap® Abrader's Main Characteristics

#### Functional Characteristics:

- Vacuum control and independent crystals flow.
- Digital indication of vacuum and crystal flow.
- Combined of crystal and filler disposables.
- Different glass tips (boron-silicate) recyclables.
- Two pieces of metallic hands sterilized.
- The smallest cost per treatment.
- No maintenance required.
- The most silent and powerful equipment in the market.

#### Technical Characteristics:

- Weight: 18 Kg.
- Input: selectionable between 220/110v.
- Frequency: 50/60 Hz.
- Power: 250 Watts.
- Vacuum level: 650 maximum mmHg.
- Caudal: 50 liters per minute.
- Crystal load capacity: 1 Kg. maximum.
- Crystal filters recover: disposables 1 Kg. maximum.

[Vendor Information](#)