



Aruba Health Distributors

Aruba Health Distributors

BC-601 Diamond

[Call for Pricing](#)

Length Width Height G.W. Power Voltage 37cm36cm 30cm 5.5kg 45W110V/220V

This unit provides a non-surgical skin refinish procedure, by using diamond peeling heads to abrade or rub off the top skin layer, then vacuuming out the particles along with any dirt and dead skin back up. This procedure removes skin debris, imperfection, blemishes, wrinkles and unwanted pigmentation on the skin.

[Vendor Information](#)