



Aruba Health Distributors

Aruba Health Distributors

BC-909 DUO Peel

[Call for Pricing](#)

Length Width Height G.W. Power Voltage 53cm 36cm 43cm 25kg 200W 110V/220V

This unit provides a non-surgical skin refinish procedure, by using sterile micro-crystals or diamond peeling heads to abrade or rub off the top skin layer, then vacuuming out the particles along with any dirt and dead skin back up. This procedure removes skin debris, imperfections, blemishes, wrinkles and unwanted pigmentation on the skin. It has a special feature, which is the ability to vacuum back up the remaining crystals on the skin.

[Vendor Information](#)